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1. **Produkt Name**

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www.hp-dent.com

IFU PMMA Blank

1. **Produkt**
 InfinityTemp/ClearTemp – PMMA for long term temporaries
2. **Composition**
 Polymethylmethacrylate and pigments.
3. **Range of dimensions of inorganic filler particles**
 Product does not contain inorganic fillers.
4. **Product Type**
 Color: Multilayer / Transparent
5. **Product specifications / Ø mm × height mm:**
 – 98,5x16; 98,5x20
 – 95x16; 95x20
6. **Warnings**
 After wearing this product, try to eat less irritating food, such as acid, hot, spicy, etc.,
7. **Contraindications**
 Avoid usage when allergic to acrylic resin.
8. **Wearing period**
 3 years
9. **Storage**
 Store in a clean and dry environment
10. **Packaging**
 Cardboard box
11. **Shade guide Recommendation**
 Vitapan classic
12. 
 Aidite (Qinhuangdao) Technology Co., Ltd.
 No.9 Dushan Road, Economic and Technological Development Zone,
 Qinhuangdao City China
13. **EC REP**
 MedNet EC-REP GmbH
 Borkstrasse 10, 48163 Münster, Germany
14. **Date of first issue of IFU**
 2020-05-09


15. **Product / Mechanical Properties**

E-Modulus	> 2,3 MPa
Flexural strength	≥ 120 MPa
Density	≥ 1,17 (g/cm ³)
Water sorption	≤ 20 (µg/mm ³)
Solubility	≤ 0,8 (µg/mm ³)

16. **Indications InfinityTemp**
 - Anterior/posterior crowns fixed on implants
 - Full-anatomic anterior/posterior crowns
 - Inlays, Onlays, Veneers
 - Bridges, full-arch restorations

- Indications ClearTemp**
- Splints
 - Surgical guides
 - Castables / Pressables (burns-out ash free)

17. **Prep-standards / Notes on design**

18. **Milling precautions**

Maintain equipment regularly (calibrate, clean, lubricate).
 Replace wear and tear parts in a timely manner.
 Before milling, check whether the milling strategy is correctly installed.

Check whether the disk holders are tightened to prevent the blank from loosening during the milling process, to avoid inaccurate results.
 Please note cutters once used for PMMA can no longer be used to mill ZrO₂.

PMMA will produce a large amount of chips, when milled dry it is suggested to use a sufficient vacuum cleaner.
 When milling InfinityTemp (multi-layer), pay attention to place the translucent side upwards.

19. **Customizing**

Using the correct tools for grinding and polishing restorations adjustment is essential. If you use an inappropriate grinding tool, grinding traces will appear, or the temperature might get too high and cause damage.

Use the following steps to finish the surface of the restoration:
 Cut the connecting rods with a cross-cutted burr.
 Polish the restoration by using the „1-step only“ P3.Polishers to create a smooth surface while preserving delicate texture and finish with SpeedPolish

While polishing please avoid to overheat the restoration. Highly concentrated alcohol cleaning solutions will cause cracks on PMMA and should therefore be avoided.
 Alternatively use mechanical cleaning methods, or the ultrasonic cleaning device.

20. **Adhesive Bonding / Connecting**

The PMMA material is easy to bond, using conventional bonding materials. When bonding to a composite follow the composite manufacturer's instructions.

21. **Contact**

HPdent GmbH
 Erwin-Dietrich-Str. 5
 D-72844 Gottmadingen
 +49 - 7731 - 3811045
 kontakt@hp-dent.com

www.hp-dent.com