**Harvard MTA – RootSeal plus EndoDirect
Bioaktiver MTA Wurzelkanal Sealer**

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**Harvard MTA – RootSeal plus EndoDirect
Bioactive MTA Root Canal Sealer**

Harvard MTA - RootSeal is a bioactive root canal sealing material based on MTA. It shows a high pH of 12. BioRootSealer can create apatite and bonds to the root canal walls. Sealing of the root canal is excellent due to virtual no shrinkage of Harvard MTA - RootSeal.

Harvard MTA - RootSeal shows an excellent radiopacity. In case of revisions, Harvard MTA - RootSeal in combination with a gutta percha Masterpoint can be removed easily.

Harvard MTA - RootSeal OptiCaps® are easily activated and the content of the capsule is easily ejected with the Harvard Applier OptiCaps®. Mixing of Harvard MTA - RootSeal OptiCaps® (mixing time 30 seconds) is achieved by a high frequency mixer with about 4,300 oscillations/minute.

1. Indications

Root canal sealer (only in combination with a gutta percha Masterpoint)

2. Contraindications

Not known

3. Side effects

Not known

4. Activation and Mixing

See "Instructions for the activation and mixing of Harvard OptiCaps®"

Mixing time for the Harvard MTA - RootSeal OptiCaps® is 30 seconds.

Attention:

Avoid lag times between the processes of activation, mixing and application as the material is in the process of setting and lag times may impair or prevent application of the material.

To prevent dehydration during setting, apply Harvard MTA - RootSeal intraoral immediately after mixing.

Working time of Harvard MTA - RootSeal is **>10 minutes** (at 23°C/74°F) and setting time in mouth is **< 60 minutes**.

5. Application:

Use of rubber dam is recommended

Prepare and clean the root canal (e.g. with Harvard Glide & Clean).

Dry the canal using paper points.

Check the Harvard EndoDirect syringe for tightness of the endo tip.

Note:

Use only the fitted flexible endo tip in combination with the Harvard EndoDirect syringe!

Mark the working length with the silicone stop. Remove the plunger from the EndoDirect syringe.

Mix Harvard MTA - RootSeal as previously described under point 4 and transfer the mixed material into the EndoDirect syringe.

Inject passively the Harvard MTA - RootSeal into the root canal to the working length. Filling should occur from bottom upwards.

Apply the gutta percha master point into the canal to the working length.

Using a heating device to cut the remaining excess of the gutta percha at the coronal part. And condense the gutta percha for a better sealing with a plugger.

Seal the cavity with a temporary filling material.

Additional remarks

- After application handle the placed MTA cement carefully.
- Store Harvard MTA - RootSeal in the sealed packaging in a dry place prior to use.
- Intraoral application of Harvard MTA - RootSeal must be done immediately after mixing to prevent dehydration during setting.
- Harvard MTA - RootSeal can cause discoloration. Use Harvard MTA - RootSeal only in the root canal.

Storage

Store Harvard MTA - RootSeal in a dry place at 2 – 25 °C (36 – 77 °F). Do not use after expiry date.

Harvard OptiCaps® and Harvard EndoDirect syringe are for single use only.

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

Limitation of Liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.











For dental use only!**Store product out of reach of children!**

Keep the instructions for use for the duration of the application.

Order-No. Article

- | | |
|---------|---|
| 7081511 | Harvard MTA - RootSeal 2 OptiCaps® ea. 0.25 g, ea. packed in an aluminium pouch, 2 EndoDirect syringes with flexible endo tip |
| 7092000 | Harvard Applier OptiCaps® |

Instructions for the activation and mixing of Harvard MTA - RootSeal OptiCaps®

	<p>Harvard OptiCaps® before activation.</p>		<p>Check the Harvard EndoDirect syringe for tightness of the flexible endo tip. Mark the working length with the silicone stop.</p>
	<p>For activation of the Harvard OptiCaps®, press the plunger on a hard and plane surface to the end into the capsule.</p>		<p>Remove the plunger from the syringe.</p>
<p>Click 1x</p> 	<p>Insert the Harvard OptiCaps® into the Harvard Applier OptiCaps® and click once. Note: The plunger must be at the same level as the bottom of the capsule.</p>		<p>Transfer the mixed material immediately into the Harvard EndoDirect syringe and apply. Unlock the Harvard Applier OptiCaps® (press button A) and remove the Harvard OptiCaps®.</p>
	<p>Insert the Harvard OptiCaps® into a mixer, close lid and mix immediately for 30 seconds at 4,300 oscillations/minute.</p>		<p>Place the plunger in the Harvard EndoDirect syringe.</p>
<p>1</p>  <p>Click 2x</p> <p>2</p> 	<p>Insert the Harvard OptiCaps® into the Harvard Applier OptiCaps®. Remove the pin from the nozzle. If you leave it, the capsule can burst. Pull the lever 2 times (2 clicks) to prime the the Harvard OptiCaps®.</p>		<p>Application of Harvard MTA - RootSeal into the root canal.</p>

Note:

Use only the fitted endo tip in combination with the Harvard EndoDirect syringe!

Only with the Harvard capsule applier (Order.-No. 7092000) the optimal amount of mixed material is guaranteed.

For the selection of a suitable capsule mixer, our sales and marketing colleagues are gladly available to you.