

## ERGO Acryl

Instrument for veneering composite




Double-sided modelling instruments ideal for acrylic facings.

### Advantages

- Three individual flexible blades, 0.2 mm thick.
- Flexible blades are secured to the instrument handle by a threaded connection and can be easily exchanged, as required.
- Ideal for applying and modelling modern facing acrylics.

### Ordering information

|   |               |
|---|---------------|
| 1 ERGO Acryl Instrument, No. 1, 1 piece | No. 1052 1100 |
| 2 ERGO Acryl Instrument, No. 2, 1 piece | No. 1052 1200 |
| 3 ERGO Acryl Instrument, No. 3, 1 piece | No. 1052 1300 |
| A Flexible blade tip, No. 1, 2 pieces   | No. 1052 1110 |
| B Flexible blade tip, No. 2, 2 pieces   | No. 1052 1210 |
| C Flexible blade tip, No. 3, 2 pieces   | No. 1052 1310 |



## ERGO Wax

Wax modeling instrument



ERGO wax instruments are characterised by the wide variety of applications for which they can be employed and by their individual adaptability. The handle design is different for each instrument, thus making them unique.

### Advantages

- Instruments can be ground according to individual needs.
- Probes, blades and scrapers made from handground stainless steel.
- Fine, anodised handles, comfortable to hold, heat-resistant and insulated.

### Ordering information

|                                      |               |
|--------------------------------------|---------------|
| 1 ERGO Wax Instrument No. 1, 1 piece | No. 1034 2001 |
| 2 ERGO Wax Instrument No. 2, 1 piece | No. 1034 2002 |
| 3 ERGO Wax Instrument No. 3, 1 piece | No. 1034 2003 |
| 4 ERGO Wax Instrument No. 4, 1 piece | No. 1034 2004 |
| 5 ERGO Wax Instrument No. 5, 1 piece | No. 1034 2005 |

