



Harvard CalciumHydroxide LC

Calcium hydroxid liner, light cure

The Harvard Calcium Hydroxide LC is suitable for indirect pulp capping and as a liner for fillings. The high pH value supports the antimicrobial effect.

Properties and advantages

- Suitable for indirect pulp capping and as liner for filling procedure
- High pH level with antimicrobial effect
- Optimal flow
- Radiopaque
- Light cure, set on demand

Indirect pulp capping.



Article information	Order no.
Harvard CalciumHydroxide LC 2 ml syringe, 3 needle tips	7071263
Harvard NeedleTips H18 refill bag with 50 needle tips	7095158



HARVARD®

Marke und Qualität seit 1892

Harvard Dental International GmbH
 Margaretenstr. 2 - 4, 15366 Hoppegarten, Germany
 Phone: + 49 (0) 30/99 28 978-0
 Fax: + 49 (0) 30/99 28 978-19
 info@harvard-dental-international.de

11024_EN_10/18

Liability is excluded for all printing errors and omissions. Before using our Harvard products, the respective directions for use should be noticed in every case. All measurements are internal measurements of Harvard Dental International.

Harvard Distribution Partner.